

Quick Test™ - Inhibitor

Quick Test™ - Inhibitor Strips

Fast and accurate testing of Inhibitor concentration in hydronic heating systems.

- Easy to use
- Fast and accurate testing
- Onsite testing with no initial set-up requirements
- Results in a matter of seconds
- Cost effective alternative to other solutions on the market
- Available on airtight waterproof containers
- Master packs of 18

Physical and chemical properties

Appearance:	Plastic Strip
Colour:	White
Odour:	N/A

Product uses

Fast and accurate testing of Inhibitor concentration in hydronic heating systems.

How to use

1. Drain small sample of system water from the system
2. Dip the Strip in water for 1 second
3. Shake off excess liquid
4. Leave for 10 seconds as the water soaks into the pad
5. After 10 seconds compare result with color scale on the side of the pot to indicate inhibitor concentration level

Environmental

ADEY is a carbon neutral company and cares about the environment. The company's award-winning product innovations contribute significantly to reduced carbon emissions, energy and water usage.

Quick Tests are classified as non-hazardous and non-toxic.

Quality

ISO 9001

Packaging, handling and storage

Available in air-tight plastic container.

Keep out of reach of children. Avoid contact with eyes and skin. In case of contact with eyes, bathe the eye with running water for 15 minutes. In case of contact with skin, wash with plenty of soap and water. Store in a cool, well-ventilated area. Dispose of in accordance with local regulations.

DIRECTIONS
Drain small sample of system water and dip test strip for one second. Shake off excess liquid. After 10 seconds, compare result with color scale below.

Optimal
Too weak

Inhibitor Quick Test™

Fast and accurate testing of MC1+® inhibitor concentration
Results within 10 seconds
Contains 50 test strips

ADEY®

Test
adey.com/us
Product ref: CH2-03-06189
M012210 LB2-01-06190

5 060106 374080